

*initial*™



*Restoration fabricated by Tim Tyndall, CDT*

# INTRODUCING THE GC INITIAL ZIRCONIA DISK

EXTREMELY HIGH QUALITY ZIRCONIUM DISKS

Use Initial IQ  
Lustre Pastes NF  
to achieve ideal  
aesthetics.



(CIP) Cold Isostatic Pressed for optimized physical properties. Available in two types for your convenience.

- **STANDARD TRANSLUCENCY**  
Ideal for veneering
- **HIGH TRANSLUCENCY**  
Ideal for full contour monolithic restorations
- Available in 12, 14, 18, & 25mm

**GC**